

Piloting & Traceability with wireless probes for continuous Chilling on racks

► The ACFRIVISU System



Piloting your wireless probes on PC

Associated to your cell, allows an identification of the products at the entry, an audiovisual signal alerts when T°C is reached at core.

The ACFRISOFT software allows a continuous recording of the T°C.

The ACFRILIST software is a smart extension of this system allowing to transfer a list of products from the PC to the ACFRIVISU.

ACFRIVISU

- Real-time display on HF probes (wireless probes)
- Identification of the products at the entry of the racks
- Recording on internal memory
- Export of datas on PC
- Edition of curves on the display
- Fits with ACFRI Soft
- Display of 1 to 8 wireless probes
- Audio signal for core T°C on HF probe



► ACFRI SOFT

- Recording & storage of the datas on PC
- Alarm level adjustable
- Automatic display of the cycle duration
- Can be connect with any kitchen equipment
- Printing & storage of the curves on image viewer

► HF PROBE (wireless probe)

- New design with clear lid
- Several types of probes (reinforced, thin,...)
- Lithium batteries, rechargeable
- Charge & connexion indicator
- Battery charger in optionoption